ShenZhen Century Health Medical Technology Co., Ltd.(SZCHMT)

Address: Room 3807B, Saige Square, HuaQiang Rd, Futian, ShenZhen, GuangDong, 518031, China Website: www.magic-genius.com Tel: +86 755 8300 3508 Mobile: +8613809890887

iCare2 9500 Particalute Respirator

Table of Contents

• Page 2:

iCare 9500 Quicky Feature

Page 3 ~4:

Product User Instructions

• Page 5~7:

China Registered Record FDA Registered Record SGS Audit ISO13485 Certificated

• Page 8 ~10:

Appendix I: GB2626 Declaration of Conformity

Appendix II: EC Declaration of Conformity

Appendix III: CE Certificate

• Page 11~17:

Appendix IV: KN95 GB2626-2006 Test Report from SMQ --- CMA & CNAS Authorized Lab

• Page 18~22:

Appendix V: KN95 GB2626-2006 Test Report from FZTTJ---CMA & CNAS Authorized Lab

ShenZhen Century Health Medical Technology Co., Ltd.(SZCHMT)

Address: Room 3807B, Saige Square, HuaQiang Rd, Futian, ShenZhen, GuangDong, 518031, China Website: www.magic-genius.com Tel: +86 755 8300 3508 Mobile: +8613809890887

iCare2 9500 -01: Big size-:11*16cm after folding, More fitting for Europen & American;

iCare2 9500-03: Big Size 11*16cm with Valve, More Freeing for Breathable

Product Grade:	KN95 GB 2626-2006, or FFP2 CE EN149:2001+A1:2009
	1st Layer Breathable nonwoven fabric ,
Material:	2nd+3rd Layer melt- blown filter,
Waterian.	4th Layer Spunlaced non woven fabric,
	5th Layer Skin-friendly non woven fabric
Certificated:	GB2626-2006, FFP2 CE EN149, FDA510K exempt
Layers:	5 layers Meltblown fabric >95%
Features:	3D fitting design, no stitching, skin friendly and breatable clothing
Use Range:	Personal protection, Professional antiviral
Packing:	10pcs*100Pack= 1000pcs/carton
Carton Size:	55x39x36cm ,Net weight 8.5kg, Volume weight 15.5kg



MAGIC GENIUS KN95 GB 2626-2006 Particulate Respirator



MAGIC GENIUS

KN95 GB2626-2006

Particulate Respirator

Foldable (For single use only)

*Picture is for reference only,actual product may have differences.





KN95 Particulate Respirator

Filtration efficiency $\geq 95\%$

Application

It is intended for respiratory protection against non-oil particles and reduce user exposure to certain airborne particles including those from sanding, sawing, grinding, sweeping, bagging, or from processing minerals, coal, iron ore, flour, and certain other substances.

Use in Schools, high-speed rail, airports, public transport and other relatively dense palces and family environment.

Warning



- •Do not modify or abuse this product.
- •This product can not produce oxygen and don't use in anoxic environment.
- •This product best continuous use time is less than 8 hours.
- •The mask does not eliminate the risk of contracting any disease or infection.
- •Improper use may lead to illness and even death.
- •Some materials in direct contact with skin may cause allergic reactions.
- •For adult use only.
- •Do not wash.

MAGIC GENIUS KN95 GB2626-2006

Particulate Respirator

Directions

- 1, Put the face mask with the nose clip on the upper side and the white face on the inner side. The rubber bands at both ends should be hung on the ears.
- 2, Gently press the nose clip to make it fit with the bridge of the nose then press the nose clip to pull the lower end of the mask to the lower jaw.
- 3, Arrange the wearing effect to the mask to cover the bridge of the nose to the lower jaw for the best protective effect.












Product composition



Specifications

Name	Parameter			
Product name	KN95 Particulate Respirator			
Model	iCare2 9500			
Size	11*16cm(±0.5m) After Folding			
Material	5Layers :3Layers Spunbond nonwovens, 2 Layers melt-blown nonwovens, elastic earloops			
Storage	Store in a cool, dry, clean place away from fire.			
	Store between 5 $^{\circ}$ C ~40 $^{\circ}$ C with relative humidity below than 50%.			
Period of Validity	1 year			
Packaging specification	10pcs/bag; 1000pcs/carton;			

Standards

Magic Genius KN95 mask complies with China standard GB2626-2006 Self-Priming, Filter-Type Mask for Respiratory Protection This standard specifies

the technical requirements, testing methods and manufacturing procedures of the mask.


ShenZhen Century Health Medical Technology Co., Ltd.(SZCHMT)

Address: Room 3807B, Saige Square, HuaQiang Rd, Futian, ShenZhen, GuangDong, 518031, China Website: www.magic-genius.com Tel: +86 755 8300 3508 Mobile: +8613809890887

China Registration No.: 914403003196161065



USA FDA Registration No.:3016690413







Report Number: QIP-ASI202466

Audit Date : 08 May.,2020

This report is issued by Focus Technology Co., Ltd. (Made-in-China.com) and the supervising inspectorate (SGS-CSTC Standards Technical Services Co., Ltd.) to confirm that:

Company Name : ShenZhen Century Health Medical Technology Co., Ltd.

深圳市世纪健康科技有限公司

Showroom : <u>Http://www.magic-genius.com</u>, <u>http://mfjl01.en.made-in-china.com</u>
Address : Room 2605, Block A, WorldTrade Plaza, Futian District, Shenzhen,

Guangdong, China

Product : Face Mask, Protective Cloth, Goggles, Face Mask Production

Machine, Eye Massager, Thermometer

has been on site audited for the Following Scope of Activity


1. General Information

- 2. Foreign Trade Capacity
- 3. Management System and Product Certification
- 4. Energy Saving and Emission Reduction
- 5. Photos

Audited Supplier

General Comments:

ShenZhen Century Health Medical Technology Co., Ltd. is an authorized company with China Type II Medical relative products operations, it was established in 2014, the Sales office is located in Room 2605, Block A, WorldTrade Plaza, Futian District, Shenzhen, Guangdong, China. The company occupies an area of more than 520 square meters. Shenzhen Century Health Medical Technology Co., Ltd. has its own brand and obtained many patent certificates. They have passed ISO9001:2015 / ISO13485:2016 certification in 2020, and obtained FDA certificate and GB2626-2006 test report for mask. ShenZhen Century Health Medical Technology Co., Ltd. has successful foreign trading experience in Europe, South America, North America, Africa, Australia and Asia.



SGS-CSTC Standards Technical Services Co., Ltd.



Unless otherwise agreed in writing, this document is issued by the Company under its General Conditions of Service accessible at http://www.sgs.com/en/Terms-and-Conditions.aspx Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon is solely limited to visual examination of the safely and readily accessible portions of the consignment and reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing /inspection report & certificate, please contact us at telephone: (86-755) 8307 1443 or email: CN.Doccheck@sgs.com Page No.: 1 of 13

| SF, B4 Bulding, Zijn (Baixia) Technology (ncutetion Special Park, No.1 Guarghua Road, Nanjing, China 210014 t (86-25) 87128385 f (86-25) 83407773 www.sgsgroup.com.cc

5所,B4Bulding, Zijin (Baixia) Technology Incubation Special Park, No.1 Guanghua Road, Nanjing, China 210014 t (86-25) 87128385 f (86-25) 83407773 www.sgsgroup.com.cn 中国・南京・秦淮区光华路1号紫金白下创业特别社区4号楼B幢5楼 邮编: 210014 t (86-25) 87128385 f (86-25) 83407773 e e as.cn@sgs.com



This is to certify that the Quality Management System of

Shenzhen Century Health Medical Technology Co., Ltd.

Business License Number: 914403003196161065 Registered Address: 3807B, SEG, Huaqiang Road, Futian District, Shenzhen, China Audit Address: Room 2605, block A, World Trade Plaza, Fuhong Road, Futian District, Shenzhen, Guangdong Province, China

applicable to

Operation on personal health protection articles, health care products and type II medical devices.

has been assessed and registered by POSI against the provisions of

ISO13485:2016

This registration is subject to the company maintaining a quality management system, to the above standard, which will be monitored by POSI.

Please consult the website: www.posicert.com

The certificate information is also available on the CNCA official website: http://cx.cnca.cn.



General Manager

han

Certificate Registration No: 381200024R0S

Initial issue date: 2020.04.24 Issue date: 2020.04.24 Valid until: 2023.04.23



ShenZhen Century Health Medical Technology Co., Ltd.(SZCHMT)

Address: Room 3807B, Saige Square, HuaQiang Rd, Futian, ShenZhen, GuangDong, 518031, China Website: www.magic-genius.com Tel: +86 755 8300 3508 Mobile: +8613809890887

Declaration of Conformity

In Accordance with the Directive Issued by

Standardization Administration of China

Respiratory Protective Equipment -- Non-Powered Air-Purifying Particle

GB 2626-2006

SHENZHEN CENTURY HEALTH MEDICAL TECHNOLOGY CO., LTD.

Hereby declares that

Part Name: Fold-able respiratory protective disposable face mask

Model number: KN95

Are in conformity with the applicable parts of the following directives:

Standardization Administration of China

Respiratory Protective Equipment -- Non-Powered Air-Purifying Particle

GB 2626-2006

In addition, I hereby declare that the part named above has been designed to comply with the relevant sections of the above referenced specifications. The unit complies with all applicable essential requirements of the directives.

For and on behalf of

深圳市世纪健康科技有限公司 SHENZHEN CENTURY HEALTH MEDICAL TECHNOLOGY CO.,LTD. Date of issue:

releates

2020.03.28 Stamp and Signature of authorized ture(s)

Declaration of Conformity

MANUFACTURE: SHENZHEN CENTURY HEALTH MEDICAL TECHNOLOGY CO.,LTD.

ADDRESS: 3807B, SAIGE SQUARE, HUAQIANG ROAD, FUTIAN, SHENZHEN,

GUANGDONG, CHINA 518031

WE DECLARE THAT THE FOLLOWING PRODUCT:

KN95 Particulate Respirator

Models No: iCare2 9500, KN95-01, KN95-02

DESCRIBED ABOVE IS IN CONFORMITY WITH THE FOLLOWING DIRECTIVE(S):

PERSONAL PROTECTIVE EQUIPMENT 2016/425

RELEVANT STANDARD(S):

EN149:2001+A1:2009

THE EU DESCRIBED ABOVE HAS BEEN TESTED BY US WITH THE LISTED STANDARDS AND FOUND IN COMPLIANCE WITH THE COUNCIL REGULATION ON PERSONAL PROTECTIVE EQUIPMENT DIRECTIVE (EU) 2016/425. IT'S POSSIBLE TO USE CE MARKING TO DEMONSTRATE THE COMPLIANCE WITH THIS PERSONAL PROTECTIVE EQUIPMENT DIRECTIVE.

For and on behalf of

深圳市世纪健康科技有限公司 SHENZHEN CENTURY HEALTH MEDICAL TECHNOLOGY CO.,LTD.

Date of issue:

Plelentas

2020.03.28

Stamp and Signature of authorized Signature(s)



Certificate

ICR Polska/VC/M6330131/A1

Name and address of Best Source International Ltd.

certificate owner: 05rm,26F, A Block,WorldTrade Plaza, Fuhong Road, Futian, Shenzhen,China

Name and address of Shenzhen Century Health Medical Technology Co.,Ltd.

manufacturer: 3807B,SaiGe Square,Huaqiang Road, FuTian, Shenzhen,GuangDong,China

Product name: Infections Prevention Face Mask (Particulate Respirator)

Product types: KN95, KN9500, FFP2-9500, KN9510, KN9520, KN9530, KN9550, KN9600, KN9610,

KN9800, KN9900, iCare2-9500, iCare2-9600, iCare2- 9700, iCare2-9800, iCare2-9900

Product trademark: Magic Genius

• This certificate confirms that the product meets the requirements of the following standard(s):

The certification has been carried out in accordance with individual rules and conditions agreed with the applicant. Evaluation has been carried out in accordance with:

Test report(s) No.: A20032480SR-01

Test conducted by: United Testing Technology (Hong Kong) Limited

 Certificate issue date:
 18, 04, 2020

 Expiration date:
 02, 04, 2025

- This certificate refers to the above mentioned product and its conformity in regards of above mentioned standard(s) was proven on test sample
- This certificate does not imply meeting all essential requirements, assessment of the series-production or any other restricted Notified Bodies conformity assessment procedure appropriate for the product
- This certificate holder shall use this certificate in connection to declaration of conformity and technical data relevant for the product the certificate was issued
- This certificate does not constitute legal base to affix marking on the product

Certificate modifications:

EN149:2001+A1:2009

Replaces cert ICR Polska/CE/V/M6330131



ICR Polska Co. Ltd.
Plac Przymierza 6
03-944 Warsaw
www.icrpolska.com
e-mail: icrpolska@icrqa.com



Director: Rafał Kalinowski

Warsaw, 18, 04, 2020

















检验报告 **TEST REPORT**



深圳市计量质量检测研究院 Shenzhen Academy of Metrology & Quality Inspection



NSTC 国家体育用品质量监督检验中心(广东)
National Sporting Goods Testing Center (Guang Dong)









Page 1 of 6 Pages

Sample Description: KN95 Particulate Respirator (Infections Prevention Face Mask)

Model/Specification/Grade: 11cm*16cm(±0.5cm) 5Layers

Applicant: Shenzhen Century Health Medical Technology Co., Ltd.

Applicant Address: 3807B, SaiGe Square, Futian shenzhen, Guangdong, 518031, CHINA

Date of Receipt: 2020-04-16

Test Period: 2020-04-16 to 2020-04-24

Test Location: Longhua Experimental Base

Shenzher	Academy o	Spection
(Stamp)	素 迷	ispection 究 <i>所</i>
Issue Da	ate 2020-	04-24

Approved by:

Signature:



No.: WT204030504 Page 2 of 6 Pages

Sample Information:

Sample Description: KN95 Particulate Respirtor (Infection Prevention Face Mask)

Trade Mark: Magic Genius

Model/Specification/Grade: $16cm*11cm (\pm 0.5cm)$ 5Layers

Serial/Batch No. of Sample: -----

Manufactured Date: 2020-04-10

Manufacturer: Shenzhen Century Health Medical Technology Co., Ltd.

Manufacturer Address: 3807B, SaiGe Square, Futian shenzhen, Guangdong, 518031, CHINA

Sample Quantity: 50pcs

Sample Description before Testing: Normal.

Client Information:

Applicant: Shenzhen Century Health Medical Technology Co., Ltd.

Applicant Address: 3807B, SaiGe Square, Futian shenzhen, Guangdong, 518031, CHINA

Applicant Telephone: +8618664568327 Applicant Post Code: -----

Test Information:

Date of Receipt: 2020-04-16

Applicant No.: 8227607

Environment Condition: $(18\sim25)$ °C $(30\sim70)$ %RH

Sampling Method: Delivered by Applicant

Judgment Basis: GB 2626–2006 Test Standard: GB 2626–2006

Test Conclusion:

Test result refer to next pages.

X 11 11 台村

Tested by: 林月明 社 自油2



No.: WT204030504 Page 3 of 6 Pages

Test Item	Requirement	Test Result	Conclusion
1. Appearance		(GB 2626-2006)	Conformity
(GB 2626-2006)	5.2Requirement	1#~2#	
		Conformity	
2. Filtration	NaC1	(GB 2626-2006)	Conformity
efficiency(%)	Non-oil aerosols	Unpretreated sample:	
(GB 2626-2006)	KN95≥95. 0	1#: 98.23	
		2#: 98.20	
		3#: 98.29	
		4#: 98.48	
		5#: 98.53	
		6#: 98.30	
		7#: 98.29	
		8#: 98.90	
		9#: 98.57	
		10#: 98.44	
		Pretreated sample:	
		1#: 95.53	
		2#: 95.34	
		3#: 95.30	
		4#: 95.24	
		5#: 95.06	
		KN-Series	
		Temperature: 23.0°C	
		Relative humidity: 35.5%	
		Aerosol chamber: NaCl	
		Concentration of aerosol chamber:	
		15mg/m^3	
		(Flow meter rate: 85L/min)	





No.: WT204030504 Page 4 of 6 Pages

Test Item	Requirement	Test Result	Conclusion
3. Resistance of	Total inhalation	(GB 2626-2006)	Conformity
inhalation and	Resistance≤350	Inhalation resistance:	
exhalation(Pa)		Unpretreated sample:	
(GB 2626-2006)		1# 2#	
		72. 7 80. 4	
		Pretreated sample:	
		1# 2#	
		81.7 91.0	
	Total exhalation	Exhalation resistance:	
	Resistance≤250	Unpretreated sample:	
		1# 2#	
		78. 7 73. 3	
		Pretreated sample:	
		1# 2#	
		90. 2 79. 4	
		(Flow: 85L/min)	
4. Dead space (%)		(GB 2626-2006)	Conformity
(GB 2626-2006)	≤ 1	Disposable facepiece:	, and the second
(,		0. 67	
		(Temperature: 23.5°C)	
		(Temperature: 20.00)	
5. Head harness		(GB 2626-2006)	Conformity
(GB 2626-2006)	Disposable facepiece:	Unpretreated sample:	, and the second
(92 2020 2000)	10N, continuous 10s	1#~2#: Pass	
	No slippage, breakage		
	2222200, 220010000	Pretreated sample:	
		1#~2#: Pass	





No.: WT204030504 Page 5 of 6 Pages

Test Item	Requirement	Test Result	Conclusion
6.Flammability		(GB 2626-2006)	Conformity
(GB 2626-2006)			
Afterflame time (s)		Disposable facepiece: Unpretreated sample: 1# 2# 0.0 0.0	
		Pretreated sample: 1# 2# 0.0 0.0 Displacement speed of Head mold nose	
		tip: 60mm/s Flame temperature at 20mm from the top of the burner: 812°C	

Note:

- 1. The sequence of temperature and humidity pretreatment is as follows:
- a) Place at (38 \pm 2.5) °C and (85 \pm 5)% relative humidity for (24 \pm 1) hours;
- b) Treatment in dry environment at (70 ± 3) °C for (24 ± 1) hours;
- c) Place at (-30 ± 3) °C for (24 ± 1) hours;

After the sample temperature is restored to room temperature for 5 hours, the test shall be carried out again.



No.: WT204030504 Page 6 of 6 Pages



END OF REPORT











检验报告 TEST REPORT







Scan is following us

Scan checks for authenticity



浙江省轻工业品质量检验研究院

Zhejiang Light Industrial Products Inspection and Research Institute

国家纺织服装产品质量监督检验中心(浙江)

National Textiles and Garment Quality Supervision Inspection Center(Zhejiang)

Zhejiang Light Industrial Products Inspection and Research Institute National Textiles and Garment Quality Supervision Inspection Center(Zhejiang)

Test report

Number:W202017474E page1/3

Name of Customer	Shenzhen Century Health Me Technology Co., Ltd.	edical	Address		B, SaiGe Square, Futian n, Guangdong, 518031, CHINA
Manufacturer	Address		B, SaiGe Square, Futian n, Guangdong, 518031, CHINA		
Sample information	Name of sample: KN95Particulate Respirator Characteristics of sample: White Trademark of sample: Specification/model: 11cm*16cm (±0.5cm, 5layers) Level KN95				
	Above-mentioned info	rmation by (Customer-	supplied	· · · · · · · · · · · · · · · · · · ·
The sent way of sample	Sent by customer		Sample q	uantity	30 pieces
Receiving Date of Sample	2020/05/11		Test Ca	tegory	Entrusted inspection
Date of Testing	2020-05-11- 2020-05-12				
RatingRequirements	GB 2626-2006				
Test Summary:					
See the att	ached page for the results.		Test	[签章 : Seal	
		Da	ate of App	proval: 2	020-05-12
Date	e of manufacture: 2020-04-28.	Win S			
Remarks	· · · · · · · · · · · · · · · · · · ·				

Approved by:



Test report


Number: W202017474E

page2/3

ITEM	STANDARD (KN95)	RESULT	RATING
1. FILTER EFFICIENCY(SALT MEDIUM) (GB 2626-2006 6.3)(%)	1個別人以品质量與 2000年		
- 政府所然(浙江)	≥95.0	99.7	PASS
2. RESPIRATORY RESISTANCE (GB 2626-2006 6.5 & 6.6)(Pa)		THE P	THE REAL PROPERTY.
-INSPIRATORY RESISTANCE	≤350	223	PASS
-EXPIRATORY RESISTANCE	≤250	191	PASS
3. HEADBAND (GB 2626-2006 6.11)			
· · · · · · · · · · · · · · · · · · ·	10N,hold on10s, No breakage or slippage	No breakage or slippage were observed	PASS
4. APPEARANCE (GB 2626-2006 6.1 & 6.2)		松	y. **
- IT	Meet the requirements of GB 2626-2006 standard	The appearance quality of 2 masks conforms to the standard	PASS
5.VISUAL FIELD BELOW MASK (GB 2890-2009)(°)	河业工参		
	≥60	77	PASS

Picture(s) of sample





DECLARATION

- 1. Our organization guarantees impartiality, independence and honesty of inspection, and is responsible for the results of inspection, keeping the samples supplied by the entrusting party confidential and at the same time protecting the ownership of the samples supplied.
- 2. The test report will be deemed invalid without signatures of the inspector/reviewer and authorized personnel, and the red special inspection stamp of our organization.
- 3. The test report will be invalid if it is altered. Copies of the report are invalid without the red special inspection stamp of our organization.
- 4. The test results shown in this report are applicable only to the samples provided by customers.
- 5. All the pages of the report are integral parts of the report. Our organization will not be responsible for any misunderstanding or other results caused by using separate page(s) of the report.
- 6. If there is any dissent of the report, the entrusting party shall notify our organization timely. For the mandatory inspection given by governmental administration departments, and dissent about the sample being tested or test results on the report should be dealt with in accordance with national regulations.

浙江省轻工业品质量检验研究院及附设的检验中心

The Affiliated Inspection Centers

浙江省轻工业品质量检验研究院

Zhejiang Light Industrial ProductsInspectionandResearchInstitute

国家纺织服装产品质量监督检验中心(浙江)

National Textiles and Garment Quality Supervision Inspection Center(Zhejiang)

地址: 浙江省杭州市江干区下沙路 300 号 6 号楼

Address: Building No.6, 300 XiaSha Road, Hangzhou, Zhejiang

国家家具产品质量监督检验中心 (浙江)

National Center for Quality Supervision Inspection of Furniture(Zhejiang)

浙江省室内安全及家具产品质量检验中心

Zhejiang Center of Quality Test for Indoor Safety and Furniture Products 地址: 浙江省杭州市余杭区良渚街道经一路 1 号良渚大学科技园 4 号楼

Address: Building 4 LiangZhu University Science and Technology Park, No.1,

浙江省轻工及五金产品质量检验中心

Zhejiang Center of Quality Test for Light Industry and Hardware Products

JingYi Rd, Yuhang District, Hangzhou, Zhejiang

浙江省体育用品质量检验中心

Zhejiang Center of Quality Test for Sports Products

地址:浙江省杭州市西湖区天目山路 222 号 3 号楼

浙江省杭州市余杭区良渚街道经一路1号良渚大学科技园3号楼

Address: No. 222 Tianmushan Rd., Hangzhou, Zhejiang

Building 3 LiangZhu University Science and Technology Park, No.1,

JingYi Rd., Hangzhou, Zhejiang

国家锁具产品质量监督检验中心(浙江)

National Center for Quality Supervision Inspection of Lock(Zhejiang)

浙江省锁具产品质量检验中心

Zhejiang Center of Quality Test for Lock Products

地址:浙江省杭州市塘苗路 24 号

Address: No. 24 Tangmiao Rd., Hangzhou, Zhejiang

联系电话: 0571-85122669 Telephone: 0571-85122669

E-Mail: 1137907994@qq.com

联系电话: 0571-89009556 Telephone: 0571-89009556

E-Mail: 2047699564@qq.com

联系电话: 0571-89001107 Telephone: 0571-89001107

E-Mail: wujtest@126.com

联系电话: 0571-85027738 Telephone: 0571-85027738

E-Mail: Locktest@126.com

网上业务受理/报告查询: http://www.zjttj.cn Online Business-Reception/Report Inquires: http://www.zjttj.cn 投诉电话: 0571-85023552 Complaint Tel: 0571-85023552